

BITOBOX THERMO- INSULATION SET

BITO
LAGERTECHNIK

Safe and reliable shipping of
temperature-sensitive products

KEEPS
TEMPERATURES
BETWEEN +2°C AND +8°C
FOR AT LEAST 12 HOURS

AVAILABLE IN 2 DIMENSIONS TO
SUIT BITOBOXES MB

EASY HANDLING




www.bito.com

KEEP YOUR ITEMS SAFE WITH OUR TEMPERATURE CONTROLLED SHIPPING SOLUTION


Thermo-insulation set for BITOBOXes MB

MAINTAINS
2-8°C
FOR AT LEAST
12 HOURS

Summer - In-box temperature range maintained with the thermo-insulation set



Winter - In-box temperature range maintained with the thermo-insulation set



Ideal for shipping temperature-sensitive products: pharmaceuticals, vaccines, etc.

Designed to minimise transport damage or product spoilage

Keeps temperatures between +2°C and +8°C for at least 12 hours

Safety-approved. Instructions for adaptation to summer and winter temperatures.*

* Example: At a constant ambient temperature of 20°C two cold packs are sufficient to keep the internal temperature between 2°C and 8°C for 8 hours.

5-litre thermo-insulation set 4327

- + Useful volume 5 litres
Inside dimensions LxWxH: 170x170x170 mm
- + 3 cool packs (white; Icecatch® Solid Isolated)
- + PS insert box (white; set is also available without insert box)
- + Neopor® inlay (grey)

To fit **BITOBOX** MBD43271
(dimensions LxWxH: 410x300x290 mm)

17-litre thermo-insulation set 6427

- + Useful volume 17 litres
Inside dimensions LxWxH: 170x170x170 mm
- + 3 cool packs (white; Icecatch® Solid Isolated)
- + PS insert box (white; set is also available without insert box)
- + Neopor® inlay (grey)

To fit **BITOBOX** MBD43271
(dimensions LxWxH: 410x300x290 mm)



3 cool packs
Icecatch® Solid Isolated
with interior coating to prevent
cold packs from freezing to
each other or to contents

Easy-to-clean insert box
resists temperatures of
up to 60°C

Insulating
Neopor® inlay

Perfectly adapted to
multi-purpose containers
- MBD43271 (400x300x290 mm)
- MBD64271 (600x400x290 mm)

BITOBOX MB

SAFE AND RELIABLE SHIPPING CONDITIONS

THERMO-INSULATION SET + BITOBOX MB

Get your sample box at a discount price! [Contact us.](#)



Reference number	Item	Number	Price € / pack	Total €
SAMPLE BOX Discount price*	5-litre thermo-insulation set + BITOBOX MB	1 pce	20.00	20.00
6-31363	5-litre thermo-insulation set			
6-11124	BITOBOX MB D43271 (small)			
6-31364	17-litre thermo-insulation set			
6-11120	BITOBOX MB D64271 (large)			

Terms of sale and payment: All prices are quoted in €, and exclude VAT. Orders with a net value of more than €200 are delivered free buyer's address within Germany, packaging is included. For orders with a net value of less than €200, a shipping and handling flat rate of €20 will be charged within Germany. Contact us for our terms of delivery to other countries. We reserve the right to make price adjustments to account for fluctuations in material and labour costs, technical modifications and printing errors. All deliveries are subject to the BITO General Terms of Sale and Delivery only. Our products and services exclusively address professional needs in trade and industry. Ownership of the goods delivered shall remain with BITO until full payment has been made. Payment is due within 30 days from the date of invoice net cash or with a discount of 2% within 8 days from the date of invoice. These terms shall be governed by German law. The place of jurisdiction agreed by the Parties for all disputes arising out of or in connection with the present contract shall be the court of Bad Kreuznach. * The shipping and handling flat rate does not apply to sample containers.

www.bito.com/katalog



BITO-LAGERTECHNIK BITTMANN GMBH

Obertor 29
D-55590 Meisenheim
Tel. +49 6753 122-0
Fax +49 6753 122-399
bitobox@bito.com
www.bito.com

A man with short brown hair and glasses, wearing a white lab coat, is looking down at a stack of blue plastic storage bins. The bins are labeled with white stickers. The top bin is labeled 'Station H2 Gefäßchirurgie' and the bottom bin is labeled 'Station G2-I Diabetisches Zentrum'. The background shows a laboratory or clinical setting with shelves and equipment.

IDEAL FOR SHIPPING TEMPER-
ATURE-SENSITIVE PRODUCTS:
**PHARMACEUTICALS, VAC-
CINES, ETC.**

BITOBOX MB

BITOBOX THERMO-INSULATION S

SAFE AND RELIABLE SHIPPING

IN MOST CASES, BINS & CONTAINERS USED IN ^{AUTOMATED} SYSTEMS ARE EXPECTED TO MEET THE FOLLOWING

REQUIREMENTS:

- + Low noise handling
- + Compatibility with all components of your intralogistics system
- + Huge storage capacity
- + Defined load capacity
- + Water drainage holes for use in sprinkler environments
- + Customizable bin subdivision

To ensure efficient product handling you want to make sure that your storage bins and boxes can be used in automated handling systems and on conveyors.



BITOBOX MB Medical in use at the University Hospital Dresden